SANITAS

SFT 22

- D** Ohrthermometer
Gebrauchsanleitung
- GB** Ear Thermometer
Operating instructions
- F** Thermomètre auriculaire
Mode d'emploi
- I** Termometro per orecchio
Istruzioni per l'uso
- RUS** Ушной термометр
Инструкция по применению

CE 0344

Serviceadresse:
 Hans Dinslage GmbH
Riedlinger Straße 28
88524 Uttenweiler GERMANY
Tel.-Nr.: (0) 7374-915766
Fax-Nr.: (0) 7374-920723
E-Mail: service@sanitas-online.de

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ENGLISH

Please read these instructions for use carefully and keep them for later use, be sure to make them accessible to other users and observe the information they contain.

1. Important notes

- This device is a sensitive electronic device. Please treat it carefully and do not expose it to any mechanical impacts.
- Do not expose the thermometer to direct sunlight.
- The device should be in the room in which the measurement is made for at least 30 minutes before use.
- The thermometer is NOT waterproof. For this reason, direct contact with water or other liquids should be avoided.
- Clean the measuring tip after each use with a soft cloth moistened with disinfectant.
- Check before each use that the lens is intact. If it is damaged, please contact the dealer or service address.
- The thermometer was constructed for practical use but cannot replace a visit to the doctor.
- This device is not intended for commercial or clinical use.
- The device is intended only for the purpose given in these directions for use.
- Repairs must be carried out only by authorised service agents. Otherwise the guarantee becomes void.
- There are people in whom different measurement results are obtained in the left and right ear. In order to record temperature changes, always measure in the same ear in the same person.
- If you have been lying on one ear for some time, the temperature is slightly raised. Wait a little while or measure in the other ear.
- As ear wax can affect the measurement, you should clean the ear before measuring if necessary.
- The ear thermometer may be used by children only under adult supervision. Measurement is usually possible over the age of 6 months. In infants under 6 months, the ear passage is still very narrow so that the temperature of the eardrum can often not be recorded and the measurement result displayed is too low.

- The measurement must not be taken in an ear affected by inflammatory diseases (e.g. discharging pus or secretion), after possible ear injuries (e.g. eardrum damage) or in the healing period after operative procedures. In all of these cases, please talk to your doctor.
- This device complies with EU Directive 93/42/EC.
- This device complies with standard EN 12470-5 Medical Thermometers, Requirements for Infrared Ear Thermometers (with maximum device).
- This unit is in line with European Standard EN60601-1-2 and is subject to particular precautions with regard to electromagnetic compatibility (EMC). Please note that portable and mobile HF communication systems may interfere with this unit. For more details, please contact customer service at the address indicated.
- If you have further questions about using our devices, please contact your dealer or Customer Service.
- Please read these instructions for use carefully and keep them for later use, be sure to make them accessible to other users and observe the information they contain.

2. What you should know about using this thermometer

Before using a forehead thermometer, an ear thermometer or a conventional rod thermometer, the following should be noted: The different thermometers listed above are suitable for measuring body temperature in different parts of the body: forehead thermometer: measurement only on the forehead, ear thermometer: measurement only in the ear, rod thermometer: rectal (in the back passage), axillary (under the arm) or oral measurement (in the mouth).

The temperature varies depending on the part of the body where the measurement is taken. The difference in healthy persons between different parts of the body can be between 0.4–1.8°F. The approximate normal temperature ranges are as follows:

- forehead temperature – measured with a forehead thermometer: 96.4°F to 99.68°F.
- Ear temperature, measured with an ear thermometer: 96.8°F to 100°F.
- Rectal temperature, measured with a conventional thermometer: 97.3°F to 100°F.
- Orally measured temperature – measured with a conventional thermometer: 96.8°F to 99.3°F.

i Note:

Temperatures measured with different thermometers cannot be compared with one another.

You should therefore tell your doctor or bear in mind if you are diagnosing yourself what thermometer you used to take your temperature and in what part of the body.

The temperature of a healthy person is also influenced by the following factors:


- The person's individual metabolism.
- Age (body temperature is higher in babies and small children and falls with increasing age. Greater temperature fluctuations occur faster and more often in children).
- Clothing.
- The outside temperature.
- The time of day (body temperature is lower in the morning and increases during the day towards evening).
- Preceding physical and, to a lesser extent, mental activity.

i Note:


Taking the temperature gives a measurement that provides information about a person's current body temperature. If you are uncertain about interpreting the results or if the result is abnormal (e.g. fever), you should consult your doctor. This also applies in the case of slight temperature changes if there are other symptoms of illness such as agitation, severe sweating, flushed skin, fast pulse rate, tendency to collapse etc.

3. Description of device

Front side

Rear side

4. Functions

This infrared thermometer measures the

- Temperature in the human ear,

The thermometer also offers additional functionalities:

- 10 memory locations for simple plotting of changes in temperature,
- Optical and acoustic fever alarm at temperatures over 38 °C,
- switchable to °C and °F.

5. Operation

This thermometer is designed for use in the ear (ear passage / eardrum) usually for persons aged over 6 months. In children under 6 months, the ear passage is still very narrow so that the temperature of the eardrum can often not be recorded and the measurement result displayed is too low.

Before taking the temperature, remove the protective cap and ensure that the sensor tip and also the ear passage are clean. If there are signs of acute inflammation (discharge of pus or secretion, pain), in the case of injuries or directly after operations on the ear, the non-affected side should be chosen. Otherwise, the measurement result can be incorrect. Use of the thermometer in different persons can be inappropriate in certain acute infectious diseases because of the possible spread of germs despite cleaning and disinfection. Discuss this with your doctor in the individual case. It must only be used without disposable covers.


There is a battery in it already so that it can be used immediately when operated for the first time. Switch on the thermometer with the „“ button. After a short self-test and a short beep the thermometer is ready for measuring the temperature in the ear.

5.1 Batteries

The Type 3 V CR2032 battery is already fitted and depending on usage has a life of around 3000 measurements. When the battery gets weak, the battery warning symbol appears . It is still possible to measure the temperature. Battery must be replaced. If the thermometer displays the battery symbol  and "Lo", the battery must be replaced.

1. Remove screw from the battery cover.
2. Open the battery cover and remove the original battery.
3. Take a new CR-2032 lithium battery and align it with the positive terminal (anode) facing upwards.
4. Place the battery into the battery compartment.
5. Replace the battery cover over the compartment.
6. Screw the battery cover back tightly

Used batteries should not go into domestic refuse. You are legally obliged to dispose of the batteries correctly. Dispose of them through your electrical dealer or local recycling centre.

Note: you will find these symbols on batteries containing toxic substances: Pb = battery contains lead, Cd = battery contains cadmium, Hg = battery contains mercury. The battery in this device does not contain toxic substances.

6. During use

Always ensure that the sensor is clean and undamaged. Switch on the thermometer with the „☺“ button.

6.1 Measuring in the ear

After a short self-test and a short beep the thermometer is ready for measuring the temperature in the ear. Because the ear canal is slightly curved, pull the ear back and up slightly before inserting the sensor tip. This is important to allow the sensor tip to be pointed directly at the ear drum. Insert the sensor tip carefully and press the button „☺“ for about 1 sec.

The end of the measuring time is signalled by a short beep. Now you can read the measured temperature. The smiley is also illuminated according to the test result:

based on: 😊 indicates that the body temperature is in the normal range, ☹️ is displayed at a temperature of or higher than 38.0 °C, i.e. fever alarm.

The time between measurements should be at least 10 seconds. This ensures that the measurement is more accurate.

6.2 Saving the temperatures measured

The temperature last measured, is saved automatically. There are 10 memory locations available for this.

The last temperature measurements can be recalled by pressing the „M“ button on the bottom of the thermometer.

6.3 Changing the temperature unit

You can display the temperature in degrees Celsius (°C) or degrees Fahrenheit (°F). Press the button „“ to change the measuring unit. Without releasing it press the „M“ button to switch to the desired measuring unit. All saved values are displayed in the new measuring unit.

7. Trouble shooting

Error message	Problem	Solution
	Ambient temperature below 60,8°F or above 95 °F	The ambient temperature must be between 60,8°F and 95 °F
	The temperature measured is higher than 109,4 °F	Operate the thermometer only between the specified temperature ranges. If necessary, clean the sensor tip. In the event of a repeated error message, contact your dealer or customer service.
	the temperature measured is lower than 93,2 °F	Operate the thermometer only between the specified temperature ranges. If necessary, clean the sensor tip. In the event of a repeated error message, contact your dealer or customer service.

8. Cleaning, storage, disposal

Clean the sensor tip after every use. Use a clean cloth or cotton bud that can be moistened with disinfectant, alcohol or warm water.

To clean the entire device, use a soft cloth slightly moistened with a mild soapy solution. Water must never get into the device. If water gets into the device, remove the battery immediately. Only use the device again when it is completely dry. Do not use any harsh cleaning agents.

Never immerse the device in water.

Always store the thermometer with the cap on in order to protect the sensor.

The device must not be stored or used at too high or low a temperature or humidity (see technical specifications), in sunlight, in association with an electrical current or in dusty locations. Otherwise, inaccuracies can occur.

If prolonged storage is intended, you should remove the battery.

The device must be tested for accuracy (calibrated) after 2 years. Please send the complete device to the dealer's or service address.

Please dispose of the device in accordance with the directive 2002/96/EG – WEEE (Waste Electrical and Electronic Equipment). If you have any queries, please refer to the local authorities responsible for waste disposal.

9. Technical data

Note: When the device is used outside the specifications, perfect functioning cannot be guaranteed.

We reserve the right to make technical changes to improve and develop the product.

Model	SFT 22
Measurement range	93.2 °F – 109.4 °F
Laboratory measurement precision	±0.4 °F between 95.9 °F – 107.6 °F ±0.5 °F < 95.9 °F > 107.6 °F
Clinical repeat precision	0.56 °F (children, 1 to 5 years), 0.50 °F (adults)
Measurement duration	1 to 2 seconds
Interval between 2 measurements	At least 10 seconds
Measurement units	°Celsius (°C) or °Fahrenheit (°F)
Ambient operating conditions	60.8 °F – 95 °F with a relative humidity of up to 85 % (non-condensing)
Storage conditions	-13 °F – 131 °F with a relative humidity of up to 85 % (noncondensing)

Model	SFT 22
Dimensions	54 x 44 x 158 mm
Weight	65,5 g including battery
Battery	1 lithium battery (type 3V CR-2032)
Memory	For 10 measurements
Explanation of symbols	Device classification type BF  Please read the instructions for use 

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